


## The preparation of tetraamminecopper(II) sulphate-1-water

### Aims

You are going to make crystals of a complex salt, tetraamminecopper(II) sulphate-1-water, using copper(II) sulphate, concentrated ammonia and ethanol.



This is the equation for the reaction:  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O} + 4\text{NH}_3 \rightarrow \text{Cu}(\text{NH}_3)_4\text{SO}_4 \cdot \text{H}_2\text{O} + 4\text{H}_2\text{O}$

### Methods

- Prepare these items:
  - a hot water bath in a 100cm<sup>3</sup> beaker (use hot water from a kettle)
  - 6cm<sup>3</sup> of ethanol in a 50cm<sup>3</sup> beaker (use a graduated pipette)
- Weigh between 1.4g and 1.6g of copper(II) sulphate using a  $\pm 0.1$ g balance.  
Weigh a test tube using a  $\pm 0.01$ g balance and record its mass.  
Add your copper(II) sulphate to the test tube, re-weigh using the  $\pm 0.01$ g balance and record its mass.
- Add 4cm<sup>3</sup> of water using a graduated pipette and place the test tube in the water bath.  
Stir gently to dissolve the copper(II) sulphate.
- Remove the test tube of copper(II) sulphate solution from the water bath.  
**Do the next step in the fume cupboard while wearing gloves.**  
Add, with stirring, 2cm<sup>3</sup> of concentrated ammonia solution to the copper(II) sulphate solution.
- Pour the contents of the test tube into the beaker of ethanol, mix, then cool the mixture.
- Using a Buchner funnel and flask, filter the crystals.  
Wash out your test tube with **cold** ethanol and add the washings to the Buchner funnel.  
Finally, rinse the crystals with **cold** ethanol.
- Carefully scrape the crystals off the filter paper and onto a fresh piece of filter paper.  
Cover the crystals with a second piece of filter paper and pat to dry the crystals.  
  
Note that, to get the crystals completely dry, you may need to repeatedly move the crystals to dry parts of the filter paper.
- Weigh a sample bottle and record its mass using a  $\pm 0.01$ g balance.  
Carefully transfer your crystals to the sample bottle. Re-weigh the sample bottle and record its mass.

### Analysis

- Calculate the relative formula masses of  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  and  $\text{Cu}(\text{NH}_3)_4\text{SO}_4 \cdot \text{H}_2\text{O}$ .
- Calculate the amount of copper(II) sulphate you used. Use your answer to calculate the theoretical yield (mass) of tetraamminecopper(II) sulphate-1-water your reaction should have produced.
- Use your recorded masses to calculate the mass of copper(II) sulphate and the actual yield of product.
- Calculate the % yield you obtained: 
$$\% \text{ yield} = \frac{\text{actual yield of product}}{\text{theoretical yield of product}} \times 100$$

### Evaluation

Comment on your % yield, and explain any loss or gain in mass compared with the theoretical yield.

# The preparation of tetraamminecopper(II) sulphate-1-water

## Technicians' Notes

Per class:

copper(II) sulphate-5-water,  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$

3 x 1dp balances

1 x 2dp balance

2 x kettles

250cm<sup>3</sup> ethanol

50cm<sup>3</sup> concentrated ammonia

crushed ice

filter paper to fit Buchner funnels

additional filter paper for drying purposes

Buchner funnels and flasks

water pumps

Per student:

1 x sample bottle with stopper

1 x test tube

1 x 50cm<sup>3</sup> beaker

1 x 100cm<sup>3</sup> beaker

1 x 10cm<sup>3</sup> graduated pipette

1 x pipette filler

1 x spatula

1 x stirring rod

1 x pair of gloves